

# **Mock Circulatory Loop (MCL) Repeatability Testing**

## **Tool Reference**

RST Reference Number: RST24CV14.01

Date of Publication: 02/07/2024

Recommended Citation: U.S. Food and Drug Administration. (2024). *Mock Circulatory Loop to Generate Variable Adult Heart Conditions for Evaluating Mechanical Circulatory Support Devices* (RST24CV14.01). <https://cdrh-rst.fda.gov/mock-circulatory-loop-generate-variable-adult-heart-conditions-evaluating-mechanical-circulatory>

## **For more information**

[Catalog of Regulatory Science Tools to Help Assess New Medical Devices](#)

## Disclaimer

### About the Catalog of Regulatory Science Tools

The enclosed tool is part of the Catalog of Regulatory Science Tools, which provides a peer-reviewed resource for stakeholders to use where standards and qualified Medical Device Development Tools (MDDTs) do not yet exist. These tools do not replace FDA-recognized standards or MDDTs. This catalog collates a variety of regulatory science tools that the FDA's Center for Devices and Radiological Health's (CDRH) Office of Science and Engineering Labs (OSEL) developed. These tools use the most innovative science to support medical device development and patient access to safe and effective medical devices. If you are considering using a tool from this catalog in your marketing submissions, note that these tools have not been qualified as [Medical Device Development Tools](#) and the FDA has not evaluated the suitability of these tools within any specific context of use. You may [request feedback or meetings for medical device submissions](#) as part of the Q-Submission Program.

For more information about the Catalog of Regulatory Science Tools, email [RST\\_CDRH@fda.hhs.gov](mailto:RST_CDRH@fda.hhs.gov).

## **Mock Circulatory Loop (MCL) Repeatability Testing**

Please refer to instructions outlined in [D'Souza et al., JBME, 2024 \(Section 2.4\)](#) and [Contarino et al., ASAIO Journal, 2023](#) along with detailed steps outlined below for conducting the MCL Repeatability Testing.

*Note: MCL Repeatability Testing may be conducted using one of the Recommended Test Conditions (e.g., Cardiogenic Shock) in [D'Souza et al., JBME, 2024 \(Section 3.2\)](#).*

1. Set up and initialize the MCL by executing steps 1-19 in document, '[USER MANUAL: Mock Circulatory Loop \(MCL\) Setup and Testing](#)'.  
Note: During MCL Initialization, note down the volume of BAF required to initially fill the loop.
2. Record the BAF volume used for MCL Initialization. This is the initial BAF volume.
3. Reference the MCL Control Table (step 20 in document, '[USER MANUAL: Mock Circulatory Loop \(MCL\) Setup and Testing](#)') for setting input controls corresponding to the selected test condition for MCL repeatability testing.
4. Execute steps 21-26 in document, '[USER MANUAL: Mock Circulatory Loop \(MCL\) Setup and Testing](#)'.  
Note: When adding or removing BAF from the loop to achieve the target fluid levels, note down the net volume of BAF added or removed from the initial BAF volume.  
The completion of this step defines the baseline state of the MCL.
5. Record the volume of the BAF in the loop as the baseline BAF volume.
6. Execute steps 21-26 in document, '[USER MANUAL: Mock Circulatory Loop \(MCL\) Setup and Testing](#)'.  
This concludes replicate (or trial) test #1 of the 10 replicate tests.
7. Repeat steps 1-6 above nine times to complete the 10 replicate tests. It is important to ensure that the same baseline BAF volume and input control settings are used for each test.