



Summary of Verification and Validation (V&V) Activities

Tool Reference

RST Reference Number: RST24CV14.01

Date of Publication: 02/07/2024

Recommended Citation: U.S. Food and Drug Administration. (2024). *Mock Circulatory Loop to Generate Variable Adult Heart Conditions for Evaluating Mechanical Circulatory Support Devices* (RST24CV14.01). <https://cdrh-rst.fda.gov/mock-circulatory-loop-generate-variable-adult-heart-conditions-evaluating-mechanical-circulatory>

For more information

[Catalog of Regulatory Science Tools to Help Assess New Medical Devices](#)

Disclaimer

About the Catalog of Regulatory Science Tools

The enclosed tool is part of the Catalog of Regulatory Science Tools, which provides a peer-reviewed resource for stakeholders to use where standards and qualified Medical Device Development Tools (MDDTs) do not yet exist. These tools do not replace FDA-recognized standards or MDDTs. This catalog collates a variety of regulatory science tools that the FDA's Center for Devices and Radiological Health's (CDRH) Office of Science and Engineering Labs (OSEL) developed. These tools use the most innovative science to support medical device development and patient access to safe and effective medical devices. If you are considering using a tool from this catalog in your marketing submissions, note that these tools have not been qualified as [Medical Device Development Tools](#) and the FDA has not evaluated the suitability of these tools within any specific context of use. You may [request feedback or meetings for medical device submissions](#) as part of the Q-Submission Program.

For more information about the Catalog of Regulatory Science Tools, email RST_CDRH@fda.hhs.gov.

Summary of Verification and Validation (V&V) Activities

V&V Activity	Description	Result												
Test fluid characterization	<p>Blood analog fluid (BAF), used as the test fluid in the MCL, was prepared to mimic the viscous behavior of blood. BAF characterization involved fluid viscosity and density measurements.</p> <p>Please refer to document, ‘Preparation of Blood Analog Fluid (BAF)’.</p>	<p>Dynamic viscosity: 3.95 mPa·s Density: 1.095 g/cm³</p> <p>Acceptable. BAF properties closely match those of blood.</p>												
Pressure transducer calibration	<p>The LV, LA, and Ao pressure transducers were calibrated using a graduated (in mmHg) cylindrical tube containing a 0-100 mmHg range of heights of fluid column.</p> <p>Please refer to document, ‘Pressure Transducer Calibration’.</p>	<p><u>Water:</u></p> <table border="1"> <thead> <tr> <th>Ao</th><th>LV</th><th>LA</th></tr> </thead> <tbody> <tr> <td>60 mV/mmHg</td><td>59 mV/mmHg</td><td>58 mV/mmHg</td></tr> </tbody> </table> <p><u>BAF:</u></p> <table border="1"> <thead> <tr> <th>Ao</th><th>LV</th><th>LA</th></tr> </thead> <tbody> <tr> <td>60.1 mV/mmHg</td><td>59.6 mV/mmHg</td><td>57 mV/mmHg</td></tr> </tbody> </table> <p>This is a characterization only test that is meant to verify the precision and reproducibility of the pressure sensors.</p>	Ao	LV	LA	60 mV/mmHg	59 mV/mmHg	58 mV/mmHg	Ao	LV	LA	60.1 mV/mmHg	59.6 mV/mmHg	57 mV/mmHg
Ao	LV	LA												
60 mV/mmHg	59 mV/mmHg	58 mV/mmHg												
Ao	LV	LA												
60.1 mV/mmHg	59.6 mV/mmHg	57 mV/mmHg												
Flow sensor calibration	<p>The CO (EP688) and Ao (ME-25PXN) flow sensors were calibrated using a traditional timed volumetric fluid collection method.</p> <p>Please refer to document, ‘Flow Sensor Calibration’.</p>	<p><u>Water:</u></p> <p>EP688 – N/A (<i>*can only be used with BAF</i>) ME-25PXN – 50 mV/L/min (sensor gain: 85%)</p> <p><u>BAF:</u></p> <p>EP688 – 100 mV/L/min ME-25PXN – 50 mV/L/min (sensor gain: 98%)</p> <p>This is a characterization only test that is meant to verify the precision and reproducibility of the flow sensors.</p>												

Frequency response filter characterization	<p>The filter frequency response on the pressure amplifier and flowmeter modules were characterized using a simplified flow loop.</p> <p>Please refer to document, 'Filter Frequency Response Characterization'.</p>	<p><u>Recommended filter for pressure:</u></p> <p>Pressure amplifier – 30 Hz low-pass filter (for LV, LA, and Ao pressure transducer channels)</p> <p><u>Recommended filter for flow:</u></p> <p>T402 flowmeter – 10 Hz filter (for ME-25PXN flow sensor)</p> <p>FM501 flowmeter – 30 Hz pulsatile frequency response (for EP688 flow probe)</p> <p>This is a characterization only test that is meant to ensure measurement accuracy.</p>
MCL simulations of the recommended target test conditions	<p>The MCL was used to simulate the pathophysiologic hemodynamics corresponding to the pre-defined five target test conditions. Qualitative pressure and flow pulses and quantitative cardiac indices were computed and compared against the target values for validation purposes.</p> <p>Please refer to documents: 'USER MANUAL: Mock Circulatory Loop (MCL) Setup and Testing', D'Souza et al., JBME, 2024</p>	<p>Detailed test results are reported and discussed in D'Souza et al., JBME, 2024.</p> <p>The absolute difference between the target and simulated hemodynamics range between 3 – 9 mmHg for the SAP, 2 – 12 mmHg for the DAP, 0.3 – 1 L/min for the mean CO, and 0.001 – 0.01 s for the cardiac cycle time.</p> <p>Acceptable. The simulated hemodynamics closely matched the set target values.</p>
MCL repeatability testing	<p>The repeatability of the MCL test system was characterized by conducting ten replicate (or trial) tests at one of the recommended test conditions, cardiogenic shock, and by keeping all MCL inputs constant among the ten replicate tests.</p> <p>Please refer to documents: 'Mock Circulatory Loop (MCL) Repeatability</p>	<p>Detailed test results are reported and discussed in D'Souza et al., JBME, 2024 (sub-section: MCL Repeatability Testing) and Contarino et al., ASAIO Journal, 2023.</p> <p>The maximum uncertainty (95% CI) in the pressure and CO pulses are: Ao pressure = 1.4 mmHg, LV pressure =</p>

	Testing', D'Souza et al., JBME, 2024, Contarino et al., ASAIO Journal, 2023.	10.2 mmHg, LA pressure = 11.57 mmHg, and CO = 1.17 L/min. Acceptable. Low uncertainty values demonstrate repeatability of the test system.
--	---	---